2019-nCoV Saliva Ag EASY TEST

(Immunochromatography)

Catalogue Number:

0709D6G001 0709D6G002 0709D6G005 0709D6G010 0709D6G020 0709D6G025

INTENDED USE

The Midstream is a lateral flow immunoassay intended for the in vitro qualitative detection of N-protein antigen from 2019-nCoV in human saliva specimens. The Midstream can be used for individuals with or without symptoms or other epidemiological reasons to suspect COVID-19 infection.

A positive result indicates 2019-nCoV infection, Please quarantine yourself and contact a doctor. Additional testing is necessary. Positive results do not rule out bacterial infection or co-infection with other viruses.

A negative result should be treated as presumptive. It do not rule out 2019-nCoV infection. Negative results should be considered in the context of a patient's recent exposures, history, and the presence of clinical signs and symptoms consistent with COVID-19. Contact a doctor and Confirm with a PCR test, if necessary.

For in vitro use only. Suitable for self-testing use.

SUMMARY AND EXPLANATION

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

PRINCIPLE OF THE TEST

This test uses double-antibody sandwich to legally detect the antigen of novel coronavirus (2019-nCoV) in saliva samples. During detection, the gold labeled anti-2019-nCoV monoclonal antibody in the labeling pad binds to the 2019-nCoV antigen in the sample to form a complex, and the reaction complex moves forward along the nitrocellulose membrane under the action of chromatography, it is captured by the anti-2019-nCoV monoclonal antibody pre-coated by the detection zone (T) on the nitrocellulose membrane, and finally a red color reaction line is formed in the T zone. If the sample does not contain 2019-nCoV antigen, a red color reaction line cannot be formed in the T zone. Regardless of whether the sample to be tested contains 2019-nCoV antigen, a red reaction line will always form in the quality control area (C).

MATERIALS AND COMPONENTS

Materials provided with the test kits

REF Component	0709D6 G001	0709D6 G002	0709D6 G005	0709D6 G010	0709D6 G020	0709D6 G025
Midstream	1	2	5	10	20	25
Instructions for use	1	1	1	1	1	1
Quick Reference Instructions	NA	NA	1	1	1	1

Materials required but not provided

1. Timer

STORAGE AND STABILITY

1. Store the test as packaged between 2-30°C.

2. The Test stable until the expiration date printed on the outer packing, the product will be expired after 24 months.

3. Do not use beyond the expiration date.

4. Do not freeze any contents of the test.

5. The test must remain in the sealed pouch until use.

TEST PROCEDURE

Before test, please read the instructions carefully.

- 1. Take the midstream to equilibrate to room temperature.
- 2. Open the aluminum foil bag, take out the midstream.
- 3. Insert the absorbent tip into the mouth. Make sure midstream is horizontally placed
- 4. Swab the absorbent tip in the mouth and tongue to collect oral fluid.
- Take the absorbent tip out from the mouth when the purple color move across the result window in the center of the midstream.

6. Wait for 10 minutes and read the results.

NOTE:

*When sampling, gently hold it in mouth and let saliva naturally adsorb on the absorbent tip.

*Do not eat, drink, or smoke prior to the test for at least 30 Minutes.

*Any saliva specimen is appropriate for testing but the saliva specimen collected in the morning, before mouth rinsed, eating or drinking, is recommended.



INTERPRETATION OF TEST RESULTS

This product can only perform qualitative analysis on the detection object

Positive Result:

If both C and T lines are visible within 10 minutes, the test result is positive and valid.

Negative Result:

If test area (T line) has no color and the control area displays a colored line, the result is negative and valid

Invalid Result:

The test result is invalid if a colored line does not form in the control region. The sample must be re-tested, using a new test.



INTERNAL CONTROL

The test contains a built-in internal control in the midstream. A color band appearing in the control region (C) is designed as an internal control. The appearance of the control line confirms that sufficient flow has occurred, and that the midstreamis working normally. If the control line does not appear within 10 minutes, it is considered an error in the test result and it is recommended to test again with the same sample and a new device.

LIMITATIONS

1. The result of the test should not be taken as a confirmed diagnosis, for clinical reference only. Judgement should be made along with RT-PCR results, clinical symptoms, epidemiological information, and further clinical data.

2. The Test performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture results performed on the same sample.

3. The test must be equilibrated to room temperature $(18^{\circ}C\sim26^{\circ}C)$ before used, otherwise the results may be incorrect.

4. A negative test result may occur if the level of antigen in a sample is below the detection limit of the test.

5. Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.

6. React less than 10 minutes may lead a false negative result; React more than 10 minutes may lead a false positive result.

- 7. Positive test results do not rule out co-infections with other pathogens.
- Negative test results are not intended to rule in other viral or bacterial infections.
- 9. Negative results should be treated as presumptive and confirmed with a molecular assay.
- 10. Clinical performance was evaluated with fresh samples.
- 11. Users should test specimens as quickly as possible after specimen collection.



1.Clinical Verification

The performance of Test was established with 232 sample collected from symptomatic patients, who with symptoms onset within 7 days.

	2019-nCoV Saliva Ag	Comparative RT-PCR Test Reslt			
	EASY TEST	Positive	Negative	Total	
	(Immunochromatography)	(+)	(-)	TOLAI	
	Detected Positive	108	1	109	
	Detected Negative	7	116	123	
	Total	115	117	232	
	Sensitivity	93.91%, 95% CI (87.97,97.02)			
	Specificity	99.15%, 95% CI (95.32, 99.85)			
	Accuracy	96.55%, 95% CI (93.34, 98.24)			

Positive results broken down by days since symptom onset:

Days since symptom onset	RT-PCR Positive (+)	2019-nCoV Saliva Ag EASY TEST (Immunochroma- -tography)	PPA
1	13	13	100%
2	32	32	100%
3	52	51	98.08%
4	69	67	97.10%
5	86	83	96.51%
6	102	97	96.00%
7	115	108	93.91%

Positive results broken down by CT value:

2019-nCoV Saliva Ag EASY TEST	Comparative RT-PCR Method (Positive by Ct Value)		
(Immunochromatography)	Positive (Ct<=25)	Positive (Ct>25)	
Detected Positive	69	39	
Total	70	45	
Positive agreement	98.57%	86.67%	

2.Limit of Detection

The experimental results show that for the virus culture concentration above 100 TCID_{so}/mL, the positive rate of detection is greater than or equal to 95%. For the virus culture concentration of 50 TCID_{so}/mL and below, the positive rate of detection is lower than 95%. So, the limit of detection of the Test is 100 TCID_{so}/mL.

3.Cross-reactivity

Cross-reactivity of the Test was evaluated. The results showed no cross reactivity with the following specimen.

No.	Specimen type	Conc.
1	HCoV-HKU1	10 ⁵ TCID ₅₀ /mL
2	Staphylococcus aureus	106 CFU / mL
3	Streptococcus pyogenes	106 CFU / mL
4	Measles virus	10 ⁵ TCID ₅₀ /mL
5	Paramyxovirus parotitis	10 ⁵ TCID ₅₀ /mL
6	Adenovirus 3	10 ⁵ TCID ₅₀ /mL
7	Mycoplasma pneumoniae	10 ⁶ CFU / mL
8	Parainfluenza virus 2	10 ⁵ TCID ₅₀ /mL
9	Human Metapneumovirus (hMPV)	10 ⁵ TCID ₅₀ /mL
10	Human coronavirus OC43	107 TCID50/mL
11	Human coronavirus 229E	107 TCID ₅₀ /mL
12	Human coronavirus NL63	107 TCID ₅₀ /mL
13	MERS-Coronavirus EMC/2012	107 TCID ₅₀ /mL
14	Bordetella parapertussia	10 ⁶ CFU / mL
15	Influenza B (Victoria strain)	10 ⁵ TCID ₅₀ /mL
16	Influenza B (Y strain)	10 ⁵ TCID ₅₀ /mL
17	Influenza A (H1N1 2009)	10 ⁵ TCID ₅₀ /mL
18	Influenza A (H3N2)	10 ⁵ TCID ₅₀ /mL
19	Avian influenza virus (H7N9)	10 ⁵ TCID ₅₀ /mL
20	Avian influenza virus (H5N1)	10 ⁵ TCID ₅₀ /mL
21	Epstein-Barr virus	10 ⁵ TCID ₅₀ /mL
22	Enterovirus CA16	10 ⁵ TCID ₅₀ /mL
23	Rhinovirus	10 ⁵ TCID ₅₀ /mL
24	Respiratory syncytial virus	10 ⁵ TCID ₅₀ /mL
25	Streptococcus pneumoniae	10 ⁶ CFU / mL
26	Candida albicans	10 ⁶ CFU / mL
27	Chlamydia pneumoniae	10 ⁶ CFU / mL
28	Bordetella pertussis	10 ⁶ CFU / mL
29	Pneumocystis jirovecii	10 ⁶ CFU / mL
30	Mycobacterium tuberculosis	10 ⁶ CFU / mL
31	Legionella pneumophila	10 ⁶ CFU / mL

4. Interference Substances

The test results do not be interfered with the substance at the following concentration:

No.	Interference substances	Conc.
1	Whole Blood	4%
2	Ibuprofen	1mg / mL
3	Tetracycline	3µg / mL
4	Chloramphenicol	3µg / mL
5	Erythromycin	3µg / mL
6	Tobramycin	5%
7	Throat spray (Menthol)	15%
8	Mupirocin	10mg/mL
9	Throat lozenge (Menthol)	1.5mg/mL
10	Tamiflu (Oseltamivir)	5mg/mL
11	Naphthoxoline hydrochloride nasal drops	15%
12	Mucin	0.50%
13	Fisherman's Friend	1.5mg/mL
14	Compound Benzocain Gel	1.5mg/mL
15	Cromoglycate	15%
16	Sinex (Phenylephrine Hydrochloride)	15%
17	Afrin (Oxymetazoline)	15%
18	Fluticasone propionate spray	15%

5. Precision

1. Test 10 replicates of negative and positive by using the reference materials of enterprises. The negative agreement and the positive agreement were 100%.

Test three different lots kits including positive and negative reference materials of enterprises. The negative results and the positive results were 100%.

6. Hook Effect

The Test was tested up to 1.6 \times 10 $^{\rm 5}$ TCID_{\rm so}/ml of heat-inactivated 2019-nCoV strain and no high-dose effect was observed.

PRECAUTIONS

1. For in vitro diagnostic use.

Use appropriate precautions in the collection, handling, storage, and disposal of patient samples and used test contents.

3. Use of Nitrile, Latex (or equivalent) gloves is recommended when handling patient samples.

4. Do not reuse the used midstream or saliva swab.

5. Should never open the foil pouch of the midstream exposing it to the ambient environment until the midstream is ready for immediate use.

6. Discard and do not use any damaged or dropped midstream or material.

7. Inadequate or inappropriate sample collection, storage, and transport may yield false test results.

8. Sample collection and handling procedures require specific training and guidance.

9. To obtain accurate results, do not use visually bloody or overly viscous samples.

10. To obtain accurate results, an opened and exposed midstream should not be used.

11. Testing should be performed in an area with adequate ventilation.

 Wear suitable protective clothing, gloves, and eye/face protection when handling the contents of this test.
Wash hands thoroughly after handling.

KEY TO SYMBOLS USED					
i	Consult instructions for use	~~~	Date of manufacturer		
2°C	Store at 2°C~30°C	2	Do not reuse		
Σ	Use-by date	淤	Keep away from Sunlight		
	Manufacturer	Σ	Contains sufficient for <n> tests</n>		
LOT	Batch code	REF	Catalogue number		
IVD	<i>In vitro diagnostic</i> medical device	Ť	Keep dry		
EC REP	Authorized representative in the European Community				



Room 218, Building 2, No.68, Nanxiang Road, Science City, Huangpu District, 510000, Guangzhou P.R.China Tel.: +86-020-82557192 service@dochekbio.com www.dochekbio.com



CMC Medical Devices & Drugs S.L. C/ Horacio Lengo Nº 18, CP 29006, Málaga, Spain Tel.: +34951214054 Fax: +34952330100 Email: info@cmcmedicaldevices.com

Doc No.: DC-IN-0709D6G Ver 1.0 GB Rel.: 2021/04/01